

WASHING THE SPACER

Wash the spacer once a week.

1. Take apart the spacer as directed in the package instructions.
2. Place in a sink filled with warm water, with a drop of mild dish soap. Do not put in dishwasher.
3. Let soak for 5-15 minutes and gently shake the pieces through the soapy water to clean any residue. (DO NOT USE A CLOTH OR BRUSH INSIDE THE MAIN CHAMBER.)
4. Pull the main chamber through the soapy water, shake out any excess water and place on a dish rack or towel to dry. ***DO NOT RINSE THE INSIDE OF THE MAIN CHAMBER***. You want to leave a soapy film on the inside to get rid of the static charge so the medication does not stick to the inside.
5. Rinse the remaining pieces (i.e. the mouthpiece, and any parts that are separate from the chamber) with clean water, shake off the excess water and allow to dry thoroughly on a tea towel or drying rack. Do not use paper towel – this may cause a static charge.
6. Put the spacer back together as directed.

Replace the spacer every 12-18 months.